
BiTrac Select™

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BiTrac Select™

Full Face Mask with
Standard Elbow
22mm Female

Interchangeable Elbow Compatible

Intended Use

The BiTrac Full Face mask is an accessory for use with CPAP or Bilevel systems used in the treatment of adult Obstructive Sleep Apnea (OSA) and/or ventilatory support.

Contraindications

1. This device should not be used on uncooperative or unresponsive patients.
2. This device may not be suitable for use on patients unable to remove the mask or are taking medication which may cause vomiting.

Warnings

1. Use only with a ventilator equipped with adequate alarms and safety systems for ventilator failure. This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
2. This mask does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designated circuit.
3. This device is not suitable for providing life support ventilation.
4. Any unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, drying eyes, eye pain, eye infections or blurred vision when receiving positive airway pressure or upon awakening should be reported to the clinician.
5. To minimize risks of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
6. If skin irritation develops, this may be due to excessive pressure or skin reaction or a combination of the two. Loosen the headgear and/or try another mask or forehead pad size.
7. This CPAP mask should only be used with CPAP systems recommended by your physician or respiratory therapist. A mask should not be used unless the CPAP system is turned ON and operating properly. The vent hole or holes associated with the mask should never be blocked.

Explanation of the Warning: CPAP systems are intended to be used with the special masks with connectors which have vent holes to allow continuous flow of air out of the mask. When the CPAP machine is turned on and functioning properly, new air from the CPAP machine flushes the exhaled air out through the attached mask exhalation port. However, when the CPAP machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocation. This warning applies to most models of CPAP systems.

Use of Interchangeable Elbow Feature:

1. Remove the existing elbow from the mask by holding the adapter ring and twisting the elbow counter clockwise to disengage the bayonet lock (Refer to Figure B.)
2. Gently pull the elbow from the mask. (Refer to Figure C.)
3. Gently press the desired Interchangeable Compatible Elbow into the mask adapter ring and turn clockwise until fully seated. Do not use excessive force.

Bayonet Lock References:

Operating Specifications: 5°C to 40°C
at humidity range of 15% to 95%

Storage Specifications: -20°C to 60°C
at relative humidity up to 95% non-
condensing


Figure A.


Figure B.


Figure C.

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241 USA

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Single Patient Use

Made in USA

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.

Pulmodyne is a registered trademark of Pulmodyne, Inc. BiTrac SE Select and Bilevel Compatible are trademarks of Pulmodyne, Inc.

Intersurgical UAB
Tel: + 370 367 66 011
e-mail: info@intersurgical.lt
Amonio p. 63, LT-18176 Patruškai, Lithuania

*European Union

5864F 06-14-2022