BiTrac Select™

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

Şekil B.

Şekil C.

BiTrac Select™ MaxShield™ Mask with Standard Elbow 22mm Female Interchangeable Elbow Compatible

Intended Use

The BiTrac MaxShield mask is an accessory for use with CPAP or Bilevel systems for ventilatory support.

Contraindications

- 1) This device should not be used on uncooperative or unresponsive patients.
- 2) This device may not be suitable for use on patients unable to remove the mask or are taking medication which may cause vomiting, have nocturnal vomiting, hiatal hernia or impaired cardiac sphincter function.
- 3) This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes.

Warnings

1. Use only with a ventilator equipped with adequate alarms and safety systems for ventilator failure. This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
2. This mask does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designated circuit.
3. This device is not suitable for providing life support ventilation.
4. Any unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, drying eyes, eye pain, eye infections or blurred vision when receiving positive airway pressure or upon awakening should be reported to the clinician.
5. To minimize risks of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
6. Make sure the mask and headgear are the correct size for the patient's face. If excessive discomfort is experienced, adjust the mask and head gear. If the discomfort continues, consider an alternative mask.


Operating Specifications: 5°C to 40°C
at humidity range of 15% to 95%

Storage Specifications: -20°C to
60°C at relative humidity up to 95% non-
condensing

Use of Interchangeable Elbow Feature:

1. Remove the existing elbow from the mask by holding the adapter ring and twisting the elbow counter clockwise to disengage the bayonet lock (Refer to Figure B.)
2. Gently pull the elbow from the mask. (Refer to Figure C.)
3. Gently press the desired Interchangeable Compatible Elbow into the mask adapter ring and turn clockwise until fully seated. Do not use excessive force.


Bayonet Lock References:

Figur A.

Figur B.

Figur C.

Pulmodyne Inc.
2055 Executive Drive
Indianapolis, IN 46241 USA
www.pulmodyne.com

EC REP QNET BV
Kantstraat 19
NL-5076 NP Haaren
The Netherlands

CE
2797

Single Patient Use

Made in USA

Latex Free

*European Union

EU*

Intersurgical UAB
Tel: + 370 387 66 611
e-mail: info@intersurgical.lt
Arniūnų g. 60, LT-18170, Pabradė, Lithuania

Caution: Federal Law (USA) restricts this device to sale by or on the order of a licensed physician.