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## Intended Use

The BiTrac SE MaxShield Select™ mask is an accessory for use with CPAP or Bilevel systems for ventilatory support.

## Contraindications

- 1) This device should not be used on uncooperative or unresponsive patients.
- 2) This device may not be suitable for use on patients taking medication which may cause vomiting, have nocturnal vomiting, hiatal hernia or impaired cardiac sphincter function.
- 3) This mask may not be suitable for use on patients with the following conditions: glaucoma, recent eye surgery or dry eyes.

## Warnings

1. Use only with a ventilator equipped with adequate alarms and safety systems for ventilator failure. This mask does not include an anti-asphyxia valve that would allow a patient to breathe if the ventilator failed.
2. This mask does not have an integral leak port; therefore a separate exhalation device must be part of the circuit. Utilize the equipment manufacturer's recommended exhalation port in the designated circuit.
4. Any unusual chest discomfort, shortness of breath, stomach distension, belching, severe headache, drying eyes, eye pain, eye infections or blurred vision when receiving positive airway pressure or upon awakening should be reported to the clinician.
5. To minimize risks of vomiting during sleep, avoid eating or drinking for three (3) hours before use.
6. Make sure the mask and headgear are the correct size for the patient's face. If excessive discomfort is experienced, adjust the mask and head gear. If the discomfort continues, consider an alternative mask operating properly. The vent hole or holes associated with the mask should never be blocked.

## Achieving the Proper Fit:

1. Use the size guide to measure the size of the patients face.
2. The inside cushion should rest above the eyebrow and directly below the bottom lip. Refer to Figure 1.
3. When the correct size is achieved place the proper size mask on the patient.
4. Place mask on patients face, and pull the bonnet overhead.
5. Clip the mask into place and adjust the straps. Refer to Figure 2.
6. Continue with CPAP therapy per your protocols.

## Use of Interchangeable Elbow Feature:

1. Remove the existing elbow from the mask by holding the adapter ring and twisting the elbow counter clockwise to disengage the bayonet lock (Refer to Figure B.)
2. Gently pull the elbow from the mask. (Refer to Figure C.)
3. Gently press the desired Interchangeable Compatible Elbow into the mask adapter ring and turn clockwise until fully seated. Do not use excessive force.

## Bayonet Lock References:


Figure A.

Figure B.

Figure C.

**Operating Specifications:** 5°C to 40°C  
at humidity range of 15% to 95%

**Storage Specifications:** -20°C to  
60°C at relative humidity up to 95% non-  
condensing

## Fitting Guide References:


Figure 1

Figure 2

Figure 1

Figure 2

**Pulmonary Inc.**  
2055 Executive Drive  
Indianapolis, IN 46241 USA

[www.pulmonary.com](http://www.pulmonary.com)

Single Patient Use

Made in USA

For more information regarding our trademarks  
visit: [www.pulmonary.com/ip](http://www.pulmonary.com/ip)

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Latex Free

EN-Medical Device

Caution: Federal Law (USA) restricts this device to sale  
by or on the order of a licensed physician.